


DDS Example with OpenICE

In this example I will show how to connect data from a DDS (Data distribution service for Real-Time Systems). Be aware, that this is very alpha 😊

First of all you have to install the DDS Feature:




After that, you can download the MD PnP Software from here:


<http://sourceforge.net/projects/mdpnp/>

<https://github.com/mdpnp/mdpnp>


Start the application



Choose Domain Id 0 (or remember the chosen value) and run "Start ICE_Supervisor".




You can now "Create an ICE Device Adapter"




And choose e.g. Simulated Pulse Oximeter. Call "Start Simulated Pulse Oximeter" and call "Start PO_Simulator".

Now you will get to the main screen:



When you double click on Puls Ox (Simulated) you will see:



Now back to Odysseus:

Create a new Project and add the following script:

```
#PARSER PQL
#RUNQUERY
sampleArray := DDSSOURCE({
    source = 'SampleArray',
    qosfile='${WORKSPACEPROJECT}/ice_library.xml',
    idlfile='${WORKSPACEPROJECT}/ice.idl',
    topic='SampleArray',
    topicType='ice::SampleArray',
    qoslibrary='ice_library',
    qosprofile='waveform_data',

    domain=0

})
```

There are two important files:

USER_QOS_PROFILES.xml

ice.idl

They contain information necessary for DDS. You can download the files from the following locations:

<https://github.com/mdpnp/mdpnp/blob/master/data-types/x73-idl/src/main/idl/ice/ice.idl>

https://github.com/mdpnp/mdpnp/blob/master/data-types/x73-idl-rti-dds/USER_QOS_PROFILES.xml

Copy them to the workspace.

If you changed the domain, please update domain value in the options above.

Now run the query, and show as table. You should see something like:

unique_device_i...	metric_id	vendor_metric_id	instance_id	unit_id	frequency	values	device_time	presentation_ti...	start	end
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003386j0	1466003386j0	1466003386004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003385j0	1466003385j0	1466003385005	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003384j0	1466003384j0	1466003384004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003383j0	1466003383j0	1466003383004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003382j0	1466003382j0	1466003382005	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003381j0	1466003381j0	1466003381004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003380j0	1466003380j0	1466003380004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003379j0	1466003379j0	1466003379005	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003378j0	1466003378j0	1466003378004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003377j0	1466003377j0	1466003377005	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003376j0	1466003376j0	1466003376004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003375j0	1466003375j0	1466003375004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003374j0	1466003374j0	1466003374004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003373j0	1466003373j0	1466003373005	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003372j0	1466003372j0	1466003372004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003371j0	1466003371j0	1466003371004	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003370j0	1466003370j0	1466003370745	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003369j0	1466003369j0	1466003370745	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003368j0	1466003368j0	1466003370745	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003367j0	1466003367j0	1466003370745	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003366j0	1466003366j0	1466003370745	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003365j0	1466003365j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003364j0	1466003364j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003363j0	1466003363j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003362j0	1466003362j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003361j0	1466003361j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[33648.113, 33...	1466003360j0	1466003360j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[30352.27, 302...	1466003359j0	1466003359j0	1466003370744	<null>
CD5sR9KKcpI7...	MDC_PULS_OXI...		0	MDC_DIM_DIM...	100	[37650.504, 36...	1466003358j0	1466003358j0	1466003370744	<null>

As you can see, the values attribute contains a list, and the time attributes a tuple.

Storing to a relational database

If you want to store this data in a database this is not possible directly. You have to "flatten" the data. This can be done as follows:

```

mappedValues = MAP({
    expressions = [
        'unique_device_identifiier',
        ['elementAt(time,0)', 'seconds'],
        ['elementAt(time,1)', 'nanoseconds'],
        'values'
    ]
},
sampleArray
)
table1 = UNNEST({
    attribute = 'values'
},
mappedValues
)

table1_WithNumbers = MAP({
    expressions = [
        'unique_device_identifiier',
        'seconds',
        'nanoseconds',
        ['counter()', 'pos'],
        ['values', 'value']
    ]
},
table1
)

table2 = MAP({
    expressions = [
        'unique_device_identifiier',
        ['elementAt(time,0)', 'seconds'],
        ['elementAt(time,1)', 'nanoseconds'],
        'metric_id',
        'instance_id',
        'unit_id',
        'frequency'
    ]
},
sampleArray
)

```

table1_WithNumbers is only necessary if you want to keep the order of the measurements.

The output will be as:

unique_device_identifiier	seconds	nanoseconds	pos	value	Metadata
CIMZHWsBALPediHAqjXB9paW...	1420625250	0	2901	33970.785	1420625250003 oo
CIMZHWsBALPediHAqjXB9paW...	1420625250	0	2900	33648.113	1420625250003 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2899	33228.88	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2898	32764.953	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2897	32337.092	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2896	32039.045	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2895	31951.002	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2894	32121.809	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2893	32551.34	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2892	33194.156	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2891	33980.664	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2890	34844.254	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2889	35742.566	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2888	36651.773	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2887	37568.79	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2886	38489.113	1420625249004 oo
CIMZHWsBALPediHAqjXB9paW...	1420625249	0	2885	39390.277	1420625249004 oo

and

unique_device_identifier	seconds	nanoseconds	metric_id	instance_id	unit_id	frequency	Metadata
CIMZHWsBALPediHAqjX...	1420625252	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625252003 oo
CIMZHWsBALPediHAqjX...	1420625251	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625251002 oo
CIMZHWsBALPediHAqjX...	1420625250	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625250003 oo
CIMZHWsBALPediHAqjX...	1420625249	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625249004 oo
CIMZHWsBALPediHAqjX...	1420625248	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625248003 oo
CIMZHWsBALPediHAqjX...	1420625247	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625247002 oo
CIMZHWsBALPediHAqjX...	1420625246	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625246003 oo
CIMZHWsBALPediHAqjX...	1420625245	0	MDC_PULS_OXIM_PLETH	0	MDC_DIM_DIMLESS	100	1420625245003 oo

Storing to a document store

You can store complex objects to a document store like MongoDB. For this, you have to convert the tuple to a KeyValueType. This can be done with the TupleToKeyValueTypeOperator.

See [MongoDB Feature](#) for more information on the Mongo DB feature.

```
#RUNQUERY
kv = TUPLETOKEYVALUE(sampleArray)
mongo = MONGODBBSINK({PORT = port, HOST = 'host', COLLECTIONNAME = 'collectionname', DATABASE = 'database'}, kv)
```